USER MANUAL (2 – 13)
MOBILE CONNECT SET

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

This guide gives you an overview of the basic functions of the Mobile Connect Set.

In situations where you are dependent on acoustic signals (e. g. road traffic, operating machines etc.), the hearing systems must not be set exclusively to wireless audio streaming (=Bluetooth mode), i.e. they must not be used as a wireless headset.

OPERATING ELEMENTS

- 1) Display
- 2) Programme buttons
- 3) Volume buttons
- 4) Buttons for phone operation
 - Press  to accept a phone call.
 - Press  to end or reject a phone call.
- 5) On/off switch
 - Press  for around 5 seconds to switch the Mobile-Connect remote control on/off.
- 6) Microphone opening
- 7) Keylock
- 8) a) Readout mode
 - Shortly press the button
- b) Hearing systems on/off
 - Hold down the readout button for longer than 5 seconds to switch your hearing systems off.
 - Press a programme button to switch the hearing systems on again.
- 9) Buttons for treble control
 - It is not possible to activate the treble control for Bizz 1 hearing systems.
- 10) 2.5 mm stereo socket
- 11) Socket for charger lead
 - Connect the remote control to a socket with the mains plug.
- 12) Strap with safety catch

SYMBOLS	
	Number of active programme
	Remote control is connected to phone, but is not activated
	Bluetooth connection to phone is active
 full	Remote control battery status full
 empty	
 charging	
	Volume up/down
	Treble up/down
	Keylock is activated/deactivated
	Incoming call
	Active call
	Hearing systems are switched off
	Error
	Search for Bluetooth audio signal/ Bluetooth audio signal active/ direct audio signal active
	
READOUT SYMBOLS	
	Sound and volume control for the hearing systems
	Rechargeable battery/normal battery status of left and right hearing system

HOW TO HOLD/POSITION OF THE REMOTE CONTROL

When watching television:

When telephoning:

PHONE CALLS

Before you use the remote control for the first time for making a phone call, you must:

- Activate the Bluetooth function on your phone.
- Link phone and remote control (“pairing”).

PAIRING

1) FIND AND ACTIVATE BLUETOOTH SETTINGS ON THE PHONE.

Activate the Bluetooth function on your phone (please consult your phone’s user manual to find out where to locate this function). A Bluetooth symbol should then be visible on your phone’s display when this function is activated.

2) SEARCH FOR “DEVICE DETECTION” ON THE PHONE.

Search for the “Device detection” menu item on your phone. On some phones the menu name may differ slightly. Please read point 3) first before you activate device detection on your phone.

3) ACTIVATE THE MOBILE-CONNECT CONNECTION MODE.

To activate the connection mode of Mobile-Connect, press and hold the button  on your remote control for around 7 seconds.

The remote control changes to pairing mode. You now have 2 minutes to complete the pairing process on your phone.

4) ACTIVATE THE SEARCH FUNCTION ON YOUR PHONE.

Your phone should find a device with the name "RCU Connect". This is your remote control.

Select "RCU Connect" to pair the phone and the remote control. Your phone now asks you for a PIN.

Enter "0000" as the PIN.

(If your phone requires it, confirm the connection as "secure" + "always connect automatically" on the phone.)

If the pairing fails, please repeat the previous steps again from point 3.

PHONE CALLS

- Press  on the remote control to accept a phone call.
- Speak into the microphone on the remote control.
- Press  on the remote control to end a phone call. Your hearing systems change automatically to the hearing programme previously used.
- Press  on the remote control to reject a phone call.


MOBILE-STATION (TRANSMITTER)

- 1) On/off switch
- 2) Volume control
- 3) Socket for power supply
- 4) 2.5 mm stereo socket
- 5) LED

RECEIVING AN AUDIO SIGNAL FROM AN EXTERNAL DEVICE


- Connect the Mobile-Station to the required device e. g. your television. To do this, use the cables and adapters supplied. (See diagram)
- Change to the Mobile-Connect programme with the remote control. A permanently lit blue LED confirms that the Bluetooth audio connection is established. You can now receive the sound from your external device with your hearing systems.

USING THE DIRECT AUDIO INPUT

If you want to listen to music with an MP3 player or other external device, you can use your hearing systems as a wireless headset.

- Connect the remote control to an external device e. g. your MP3 player. To do this, use the cables and adapters supplied.
- Change to the audio programme with the remote control. The remote control displays the symbol for a direct audio signal. You can now receive the sound from your MP3 player with your hearing systems.

GENERAL SAFETY INSTRUCTIONS

- Wireless devices can cause disruptions to measuring instruments and electronic devices.
- Switch off the remote control in hospitals or on planes.
- Before using the remote control with an electronic implant or life-support system, have the electromagnetic compatibility tested. Do not use the remote control if you have a pacemaker.
- Your remote control is sensitive to extreme heat, high humidity, strong magnetic fields, X-rays and mechanical stress.
- Store hearing systems, rechargeable batteries and other accessories well out of reach of babies and small children as there is a danger of swallowing and suffocation. Consult a doctor or hospital immediately if something is swallowed.
- Do not expose your remote control to extreme temperatures or high humidity, do not leave it lying in strong sunlight and do not immerse it in water.
- The wireless functionality may be impaired by an electromagnetic source (e. g. computer monitor) in the near vicinity. Remain a suitable distance away from such sources of interference.

■ If your hearing system can be controlled by another person's remote control, please take your remote control to your hearing system specialist.

■ Do not use the remote control in places where there is a risk of explosion (e. g. mines).

■ Danger of strangulation if the strap gets caught e. g. in moving parts of machines. Therefore, please use the remote control exclusively with the safety strap (the safety catch should release if the strap gets caught).

■ Some countries impose restrictions on the use of radio transmission devices. More information can be obtained from the responsible local authorities.

SAFETY INSTRUCTIONS FOR RECHARGEABLE BATTERIES

Danger of explosion if the rechargeable battery is not handled correctly:

■ Do not short-circuit the rechargeable battery.

■ Do not pierce, deform or open the rechargeable battery.

■ Do not place the rechargeable battery in the oven or on an open fire.

■ Do not immerse the rechargeable battery in water.

■ Charge the rechargeable battery at an ambient temperature of 0 °C to 45 °C.

■ Remove the rechargeable battery from the device if you are not going to use it for a long time.

CLEANING AND CARE

- Clean the remote control with a soft cloth.
- Do not use any abrasive cleaning fluids with aggressive additives such as citric acid or ammonia.

NB

Do not immerse the remote control in water. Do not clean your remote control with alcohol or benzene.

TECHNICAL DATA

1. OPERATING CONDITIONS

Ambient temperature for charging:	0°C...+45°C
Ambient temperature for operation:	-10°C...+50°C

2. POWER PACK; MODEL UE 10W-050080 SPC

Input voltage	100 – 240 VAC
Input frequency	50 – 60 Hz
Power consumption	150 mA
Output voltage	5.0 VDC
Charge current	800 mA

3. MOBILE-CONNECT

Input voltage	5 V
Charge current	500 mA
Charging time	4 h
Transmitting power	5 mW
Transmission frequency	3.28 MHz and Bluetooth
Operation period	6 h

RECHARGEABLE BATTERY – LI-ION RECHARGEABLE BATTERY

Operating voltage	3.7 V
Capacity	1000 mAh

4. MOBILE-STATION

Input voltage	5.0 V
Bluetooth transmission	Bluetooth
Operation period	5 h

DISPOSAL INSTRUCTIONS

Within the European Union, the accessories listed correspond to the European Parliament and Council Directive 2002/96/EC of 27 January 2003 on waste electrical and electronic equipment.

Later modified by Directive 2003/108/EC (WEEE.)

Do not dispose of rechargeable batteries in household waste! When they wear out, please return your rechargeable batteries to your hearing system specialist for disposal. Please send the Mobile Connect Set back to the manufacturer for disposal (address see back page).

With the CE mark, Audio Service confirms compliance with European Directive 93/42/EEC on medical devices.

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Audio Service GmbH • Zeppelinstraße 9 • D-32051 Herford